



Welcome to W M C

Walker-McGuire & Companies, Inc., your manufacturing source for products designed for distributors of Health, Sports and Medical furnishings since 1978.

Located in the heartland of America, W M C is an excellent source for domestic customers drop ship needs. Call today to join our exciting team and learn more about our products.

Content

| Introduction | 2 |
|--|-------|
| Recovery Couches | 3-7 |
| Couch Rails | 8 |
| Couch & Base description | 9 |
| Treatment Tables | 10-13 |
| Modular Taping Center w/adjustable seating | 14-15 |
| Treatment & Exam Table description | 16 |
| Modular Cabinets | 17 |
| Exam Cabinets | 18-19 |
| Storage, Lockers, Exam Tables and Carts | 20-21 |
| Modular Cabinet description | 22 |
| Laptop Workstation / Whirlpool Table | 23 |
| Office Furniture | 24 |
| Office Furniture description | 25 |
| Space Advantage Tables | 26 |
| New Products | 27 |
| Melamine & Pull samples | 28 |
| Wood Veneer finishes | 29 |
| Wood Finish care & cleaning | 30 |
| Vinyl Color chart | 31 |
| Vinyl Cleaning guide | 32 |
| Modular Taping design guide | 33 |

Frequently asked Q & A

Yes—Products are TSCA Title VI and CARB Phase 2 compliant.

Yes-Products are latex free.

Yes—Foam products meet California TB-117-2013 requirements.

Yes—Vinyl contains agents to protect against bacterial & fungal micro-organism

Yes—Wood component vendors are FSC certified.

Thank you for considering W M C for your next project.



Recovery Couches



8002: Economy couch w/chrome legs Shown: Turquoise (Pillow included) 72"L x 24"W x 19 1/2" H

Please note: 8002 & 8004 are the only couches we manufacture that does not include a spring foundation and are not available in a base configuration



9002 Flat Pillow



8004: Economy couch w/wood legs Shown: Black (Pillow included) 72"L x 24"W x 19 1/2"H

Please note: 8002 & 8004 are the only couches we manufacture that does not include a spring foundation and are not available in a base configuration.



8102: Standard couch w/chrome legs Shown: Dusty Jade w/clear cover 9151 Clear cover sold separately 72"L x 29"W x 22"H (Pillow included)

Features spring foundation



9151 Clear Covers



8104: Standard couch w/wood legs Shown: Dove (Pillow included)


72"L x 29"W x 22"H

Features spring foundation



8106: Standard couch w/standard base (Pillow included) Shown: Bronze w/Oak base 72"L x 29"W x 25"H

Features spring foundation



9004 Wedge Pillow



8108: Standard couch w/drawer base (Pillow included) Shown: Forest w/Maple base / Alum pull 72"L x 29"W x 25"H



Couches with fixed headrest

9006 Round pillow



8302: Fixed headrest couch w/chrome legs Shown: Tea Rose 72"L x 27"W x 22"H

Features spring foundation



8304: Fixed headrest couch w/wood legs Shown: Royal Blue 72"L x 27"W x 22"H

FHR

Features spring foundation



8306: Fixed headrest couch w/standard base Shown: Crocus / Walnut 72"L x 27"W x 25"H



9008 Half Round pillow



8308: Fixed headrest couch w/drawer base Shown: Black / White / Black hooded pull 72"L x 27"W x 25"H

Features spring foundation

Features spring foundation

XL Adjustable Headrest Couch



8702: XL Adjustable headrest couch w/chrome legs

Shown: Greystone 76"L x 29"W x 22"H

Features spring foundation





8704: XL Adjustable headrest couch w/wood legs

Shown: Regimental Blue 76"L x 29"W x 22"H

Couches with adjustable headrest



8502: Adjustable headrest couch w/chrome legs Shown: Burgundy 72"L x 27"W x 22"H

AHR

9004 Wedge Pillow



8504: Adjustable headrest couch w/wood legs Shown: Blue Ridge 72"L x 27"W x 22"H

Features spring foundation



Features spring foundation



9002 Flat Pillow



8506: Adjustable headrest couch w/standard base Shown: Greystone / Cherry 72"L x 27"W x 25"H



8508: Adjustable headrest couch w/drawer base Shown: Cinnabar / Walnut / Alum pull 72"L x 27"W x 25"H

Features spring foundation

Features spring foundation

Deluxe Adjustable Headrest Couch



8712: Deluxe Adjustable headrest couch w/chrome legs Shown: Black 76"L x 29"W x 24"H



Features spring foundation



8714: Deluxe Adjustable headrest couch w/wood legs Shown: Claret 76"L x 29"W x 24"H



Youth Couches with adjustable headrest



8802: Adjustable headrest couch w/chrome legs
Shown: Teal
52"L x 27"W x 22"H

Features spring foundation



8804: Adjustable headrest couch w/wood legs Shown: Cinnabar 52"L x 27"W x 22"H

Features spring foundation



8806: Adjustable headrest couch w/standard base Shown: Black / Cherry 52"L x 27"W x 25"H

AHR

8808: Adjustable headrest couch w/drawer base
Shown: Black / White base / Hooded pull 52"L x 27"W x 25"H

Features spring foundation

Features spring foundation





We've got you Covered with, NAUGAHYDE



Recovery Couches

Youth Couches with fixed headrest

FHR



8902: Fixed headrest couch w/chrome legs Shown: Forest 52"L x 27"W x 22"H

Features spring foundation



8904: Fixed headrest couch w/wood legs
Shown: Forest w/9050 Rails (sold separately) 52"L x 27"W x 22"H

Features spring foundation



8906: Fixed headrest couch w/standard base Shown: Claret / Maple 52"L x 27"W x 25"H

Features spring foundation



8908: Fixed headrest couch w/drawer base Shown: Dutch Blue / Fog / Hooded pull 52"L x 27"W x 25"H

Features spring foundation



8910: Extension Stool

May be used as an extension to couch or in a variety of ways.

Shown: Claret / Fog base (L) Dutch Blue / Maple base (R) 20"L x 27"W x 25"H

Features spring foundation



8908: Fixed headrest couch w/drawer base

8910: Extension stool

Shown: Claret / Fog / Black hooded pull

When used as a couch extension the entire length measures, 72"L x 27"W x 25"H



W M C - Model 9050 Bed rail for couches w/wood or chrome legs





9050 Couch rails

These sturdy rails can be order and installed on your existing couches, or ordered and shipped with your new couch purchase.



This illustrates the ease of installing your new Bed Rails

Begin by placing your couch upside down on a clean surface (Table or floor) with legs attached.

Place cross bar across bottom of couch and position next to couch leg with adjustment knobs facing the center of the couch. (figure 1 & 2)

Flush the wooden bracket on the cross bar with the bottom of the couch frame and secure the bracket to the frame using the $\#8 \times 11/2$ " screws provided. (a cordless screw driver will make this task effortless, however this can be done using a common #2 Phillips screws driver.)

Once the wooden mounting brackets are secured, turn the couch right side up.


Insert the outer telescoping rail into the barrel of the head end cross brace, pull the spring loaded knob on the barrel and adjust the bed rail to the top position and release the knob to secure the rail. Repeat this process for the inner telescoping rail at the foot end cross brace.

Raise or lower the bed rail until the spring loaded knob (clicks) securely into the lower or upper adjustment holes, or until the rail stop securely rest against the barrel.

Repeat these steps to install the other bed rail onto the opposite side of the couch.

Sold in sets: Each set includes 2 pre-adjusted cross bars made specifically for your couch model and 2 telescoping rail sets to secure both sides of your couch. As shown attached to WMC models 8102, 8104 & 8904 Recovery Couches.









Recovery Couches

PRODUCT DESCRIPTION & TECHNICAL DATA

RECOVERY COUCHES

Models with headrests:

Features adjustable multi-position and/or fixed headrest with paper dispenser, innerspring/no sag foundation, 2" foam and hardwood frame. Upholstered in your choice of 25 *Naugahyde vinyl colors.

*Naugahyde is a product of Uniroyal Engineered Products, LLC. (Made in the USA)

Models without headrests:

Features same as above with the exception of our Economy model which excludes a spring foundation.

Frame: Made from 21mm 15 ply Birch Plywood using ULEF (ultra-low-emitting formaldehyde) resins, which are consistently below the Phase 2 emissions standards. TSCA Title VI and CARB Phase 2 compliant!

Springs: (13) small loop 9 gauge no/sag springs, tied and covered with *DUON.

Foam: 2" flexible A-33 polyurethane foam (with range of 27 to 33 lb. density) which meets California requirement (TB 117-2013) containing NO added flame retardants.

Fasteners: Frame is assembled using 16 gauge x 2" staples and #8 x 1 1/2" pocket screws.

Vinyl: *Naugahyde developed to meet the rigorous demand of hospital, contract and institutional use. Meets the following flame retardant certification requirements: Cal 117-2013, MVSS-302, Port Authority of NY/NJ, FAA, FAR 25.853 A, BIFMA Class 1, Boston Fire Code IX-1, Fed. Spec. CID A-A-2950, NFPA 260, IMO FTP Code, Part 8. Mildew Resistant, Sulfide Stain Resistant, Oil Resistant, Cold Crack -20 degrees with advanced *BeautyGuard protective finish.

Base Panels: Made of 3/4" Melamine with the exception of the drawer sides, which are 1/2" MDF Melamine. *Funder—Melamine Panels meet or exceed American National Standard A208.1 for 1-M-2 board, ALA and NEMA LD3-1995 and meet all requirements of applicable standard CS-236-66 for 1B2 mat-formed board. (TSCA Title VI Compliant) All materials mentioned above are edged with PVC edge banding.

Hardware: Drawer slides are a premium quality bottom mount 100lb. capacity with *DELRON* rollers tested to 50,000 cycles. There are various other glued or stapled joints.

Pulls: Drawer pulls are available in oak, maple, aluminum, satin chrome wire and black or white hooded pulls.

Shipping: Each complete unit is shipped in two separate containers. The couch is wrapped in a plastic shrink film with foam protection at the corners and shipped in a double wall box. The base is shipped on a sturdy pallet with edge and corner protection in a heavy duty cardboard container.

Melamine is a *Funder America product.—*Naugahyde is a product of Uniroyal Engineered Products, LLC.















7210: Deluxe Taping Table w/2 fixed shelves.
Shown: Forest / Maple

48"L x 24"W x 36"H



7240: Space Saver Exam Table w/footrest (Not shown)

7241: Space Saver Exam Table w/foot and backrest Shown: Black / Maple 51 to 68"L x 24"W x 30"H



7220: Taping Table w/shelf (Not shown)

7221: Same as (7220) w/H-Brace Shown: Folkstone Laminate / Cherry

48"Lx 24"W x 36"H



7250: Taping Table w/center cabinet and shelves

Shown: Forest / Maple

48"L x 24"W x 36"H



7230: Bariatric Table w/back rest (Not shown)

78"L x 30"W x 30"H

7231: Bariatric Table without backrest

Shown: Regimental Blue / Maple



7260



7260: Taping Table w/Storage Cabinet, 2 doors & shelf, Shown: Regimental blue / Maple 48"L x 24"W x 36"H

Treatment Tables





7310



7270: Treatment Table w/center shelves

Shown: Sun yellow / Maple

78"L x 24" W x 30"H

7271: Same as 7270 w/back rest.

7310: Single Taping Station

Shown: Royal Blue / Walnut

44"L x 36"W x 38"H

Group these units back to back and end to end making as may stations as required






7280

7280: Taping Cabinet w/lockable storage and end shelves.

Shown: Forest / Cherry

48"L x 24"W x 36"H



7321: Table w/dual elevators and back rest

7322: Table w/standard flat top

7323: Table w/standard flat top and backrest 78"L x 30"W x 30"H





7300



7300: Dual Taping Station

Shown: Black / Maple 69"L x 36"W x 38"H

Model 7300 is a 3 piece modular cabinet that can be disassembled allowing them to be move though confined spaces.

7360: Stretch Table w/upholstered top

Shown: Black / Oak 78"L x 24"W x 19"H

7363: Stretch Table w/upholstered top

78"L x 30"W x 19"H





7100: Deluxe Table w/three drawers

Shown: Royal Blue / Maple 78"L x 30"W x 30"H

7101: Same as (7100) w/backrest.



7100



7130: Deluxe Table w/enclosed cabinet

Shown: Blue Ridge / Maple 78"L x 30"W x 30"H

7131: Same as (7130) w/backrest



7110: Deluxe Table w/Dual leg elevators

Shown: Forest / Maple 78"L x 30"W x 30"H

7111: Same as (7110) w/backrest.



7140



7140: Deluxe Table w/enclosed cabinet and sliding doors

Shown: Black / Maple 78"L x 30"W x 30"H

7141: Same as (7140) w/backrest



7120: Deluxe Table w/2 fixed storage shelves

Shown: Dove / Maple 78"L x 30"W x 30"H

7121: Same as (7120) w/backrest



7150: Deluxe Table w/Shelf

Shown: Black / Cherry 78"L x 30"W x 30"H

7151: Same as (7150) w/backrest

Treatment Tables



7160: Standard Treatment Table w/shelf

Shown: Black / Maple 78"W x 24"W x 30"H

7161: Same as (7160) w/backrest



7190: Multi-use Mat Table

Shown: Black / Maple 84"L x 45"W x 21"H

7290: Multi-use **mini** Mat Table 78"L x 30"W x 21H



7170: Deluxe Table w/nose / face cutout and shelf

Shown: Black / Maple 78"L x 24"W x 30"H

7171: Deluxe Table 78"L x 30"W x 30"H

7172: Standard Table 78"L x 24"W x 30"H w/backrest 7173: Deluxe Table 78"L x 30"W x 30"H w/backrest



7200: Space Saver Exam Table w/adjustable footrest

Shown: Royal / Maple 68"L x 24"W x 30"H

7201: Same as (7200) w/backrest



7180: Standard Treatment Table w/H-Brace

Shown: Tomato Red / Maple 78"L x 24" W x 30"H

7181: Same as (7180) w/backrest



7202: Space Saver Treatment Table w/adjustable back and foot rest

Shown: Forest / Oak 51" to 68" L x 24"W x 30"H




Modular Taping Center with adjustable seating

Modular Taping Cabinets with adjustable seating allows the trainer to easily adjust the seating heights to 21-24-27-30-33-36-

39-42".

All components are modular which allows you to build an entire room or begin with a few stations and add additional units anytime. Cabinets feature various drawer and door combinations and options such as an 18" cabinet w/trash chute & space for trash container. Seating features: plain seating—seating w/storage—seating w/individual leg pull outs and seating with an adjustable calf lift. Customize your order using our MODULAR TAPING DESIGN SHEET (page 33) or you can down load a copy from our web site. Use the worksheet to arrange your group and draw each component to get your dimensions and see how it will look before you order. These units are extremely durable, made from materials that are easy to clean and can be re-arranged to accommodate each trainer. All interiors are white with a selection of 4 solid and 4 wood grain exteriors. You will also select 1 of 25 colors from our *Naugahyde* vinyl chart found in our catalog and on our website (www.wmc-products.com). You may also want to inquire about your *logo*.



Shown are product numbers above and below the **CT (cabinet top)** This is an example of how easy this group is to design and order. Please note: Cabinets may join each other, however seating must be attached to a cabinet on each side to allow for height adjustments. You may also select one color from this chart for cabinets and a different color for seating or create all the same color!



Adding a Table extension: You may add these tables to the right, left or both ends of each group.



7650—Table end without storage,

7651—with individual leg pull out.

7652—with adjustable calf lift.

Please note: 7651 and/or 7652 will be placed on the side facing the front of each group (Shown: Maple)



7660—Table end with storage

7661—with individual leg pull out

7662—with adjustable calf lift

Please note: 7661 and/or 7662 will be placed on the side facing the front of each group (Shown: Maple)

Modular Taping Center with adjustable seating

Hutch & back selections



7630—18" Hutch 18"W x 9"D x 30"H w/1 door and 1 shelf Shown: Oak W/Black hooded pull



7634—24" Hutch 24"W x 9"D x 30"H w/2 doors and 1 shelf Shown: Oak W/Black hooded pull





7633—Seating backs
Attaches to seating items below
20"W x 3"D x 36"H
Padded area: 18"W x 22"H
Logo: submit your artwork for quote
Shown: Oak w/Black upholstery

Cabinet selections



7620—18" Base cabinet 18"W x 24"D x 42"H W/1 door & 1 drawer Provides area for trash can below trash grommet in top. Shown: Oak w/Black hooded pull



7624—24" Base cabinet 24"W x 24"D x 42"H w/1 drawer, 2 doors & 1 shelf Shown: Oak w/Black hooded pulls



7622—18" Base cabinet 18"W x 24"D x 42"H w/1 door and 2 shelves Shown: Oak w/Black hooded pull



7626—18" Base cabinet 18"W x 24"D x 42"H w/1 drawer & 1 door w/2 shelves Shown: Oak W/Black hooded pulls

Seating selections



7640—Plain seating. 28"W X 23"D x 12 1/4"H Shown: Oak w/Black upholstery



7646—Plain seating w/dual pull outs. 28"W x 23"D x 12 1/4"H With independent pulls out to support lower extremities while taping. Shown: Oak w/Black upholstery



7642 - seating w/storage. 28"W x 23"D x 12 1/4"H Provides storage under top. Shown: Oak w/Black upholstery



7648— Plain seating w/12 position adjustable calf support for taping. 28"W x 23"D x 12 1/4"H Shown: Oak w/Black upholstery



Treatment Tables

PRODUCT DESCRIPTION & TECHNICAL DATA

TREATMENT & EXAMINATION TABLES

Finish: Available in walnut, cherry, oak or natural finish with a catalyzed water borne top coat for an extremely hard and durable finish. *A gualente* stains and top coats are GREENGUARD certified products from M.L. Campbell wood finishing systems which produce superior resistance to marring, scratching, water, alcohol (200 proof) and household chemicals.

Vinyl: *Naugahyde developed to meet the rigorous demand of hospital, contract and institutional use. Meets the following flame retardant certification requirements: Cal 117-2013, MVSS-302 Port Authority of NY/NJ, FAA, FAR 25.853A, BIFMA Class 1, Boston Fire Code IX-1, Fed. Spec. CID A-A 2950, NFPA 260, IMO FTP Code, Part 8. Mildew Resistant, Sulfide Stain Resistance, Oil Resistant, Cold Crack -20 degrees with advanced *BeautyGuard protective finish.

Foam: 2" flexible D-50 polyurethane foam (with range of 48 to 50 lb. density) which meets California requirement (TB 117-2013) containing NO added flame retardants.

Frame: Made of 3/4" veneer core white Maple hardwood plywood, edged with Maple wood edge banding. Shelf made of same material. All materials are CARB Phase 2 compliant.

Legs: Made of select hardwood furniture squares, most commonly Maple or Ash. Finished dimension 2-1/4" x 2-1/4" x required length.

H-Brace: Made of select #1 or FAS kiln dried Ash hardwood. Vendors are FSC certified.

Hardware: Leg brackets are stamp zinc plated steel 3^x x 5-3/4> attached to frame with 1/4-20 carriage bolts x 4 ea. Legs attach to bracket with 3/8" x 3" lag bolt. Top is attached to frame using #8 x 1 1/2" Pocket screws, same applies to shelf. H-Brace attached using #8 x 2" screws.

Motion Models: (Models with adjustable foot and backrest.) Footrest models are fabricated using (1) pair each 3.5 x 3.5 butt hinges with a (12) position adjustable mechanism (s). Backrest models have (1) pair each 2" x 5" 10 knuckle heavy duty modified lift hinge and adjustable mechanism (s).

Cabinet Models: (Deluxe Space-Saver Exam Tables with Medical Cabinet) Euro-style design with continuous wood, aluminum or wire pulls. Door hinges are Blum® series snap-on with 110 degree opening and independent side, depth and height adjustments.

Paper Dispenser: All tables include paper, dispenser and paper cutter.

Shipping: Units will be shipped set up and packaged on a sturdy pallet (no assembly required)

*Naugahyde is a product of Uniroyal Engineered Products, LLC.











Modular Cabine

Wall Cabinets - 30"H x 12"D



Typical wall & base configuration w/Melamine laminate finish

Wall cabinet left to right Model 6010 & 6011

Base cabinet left to right Model: 6021—6031 & 6041 w/wire pulls



Typical base cabinet configuration with models 6040, 6031, and 6020 In walnut finish w/wire pull

Walnut wood finish



CUSTOM

Typical base configuration w/aluminum pulls. Left to right 6021/6041/6032 Custom maroon laminate finish, additional colors available in Formica & Wilson Art laminates



18" wide



6001: Clear Acrylic Framed Doors 24" wide



6002: Clear Acrylic Framed Doors 30" wide



6010: Solid Door 18 " wide



6011: Solid Doors 24" wide



6012: Solid Doors 30" wide

Add 1 1/2" to a group of Wall cabinets to obtain total width of finished cabinet group. Add 2" to a group of Base cabinets to obtain total width of finished cabinet group.

Base Cabinet — 36"H x 19"D



6020: Cab. w/5 Drawer 18" wide



6021: Cab. w/5 Drawer 24" wide



6022: Cab. w/5 Drawers 30" wide



6030: Cab. w/1 drawer and 1 Door, w/1 Adj. Shelf. 18" wide



6031: Cab. w/2 drawers and 2 Doors, w/1 Adj. Shelf. 24" wide



6032: Cab. w/3 Drawers and 2 Doors. 30" wide



6040: Cab. w/1 Door and 2 Adj. Shelves 18" wide



6041: Cab. w/2 Doors and 2 Adj. Shelves 24" wide



6042: Cab. w/2 Doors and 2 Adj. Shelves 30" wide



6050: Cab. with wheelchair accessible Sink. 36" wide



6051: Cab. w/2Doors and sink. 24" wide



6052: Cab. w/2 Doors and sink. 30" wide



6053: Cab. w/4 Doors and sink. 24" wide



6054: Cab. w/4 Doors and sink. 30 wide



6062: Corner base cabinet w/2 doors

36" L & R of corner

6031



6060: Corner wall cabinet w/solid doors & 4 shelves



6061: Corner wall cabinet w/acrylic doors & 4 shelves



6062: Corner Base Cabinet w/Lazy Susan & 2 shelves



MAPLE WALNUT ALMOND OAK CHERRY FOG

6021

WHITE

BLACK FOLKSTONE



Pre-configured cabinet groups allow reduction in cost by using only the components required and the options you choose. Features can be found on our Modular Medical Cabinet product description and technical data page.

Panels, doors & drawer fronts are edged in matching color 3mm PVC with rounded profile & corner rounding.

Options: 6093—Stainless steel sink w/goose neck faucet & strainer basket, available on 24" & 30" cabinets.

6096—Work surface slide out available on 18" & 24" cabinets.

6099—Locks - you may choose to lock 1 or more doors and/or drawers. (2 door cabinets require only 1 lock)

Product numbers are given to both wall & base units, allowing you to mix or choose wall and/or base groups. All groups are available **left handed**, from left wall & **right handed** (reversed from right wall), and **centered**.

When ordering please submit—Group number (s) / right or left handed / color and option (s).



6201—Wall cabinet 19 1/2"w x 30"h x 12"d w/1 door 6301—Base cabinet 26"w x 36"h x 18"d 1 drawer, 2 doors and 1 shelf. (**Left hand**)




6202—Wall cabinet 19 1/2"w x 30"h x 12"d w/1door 6302—Base cabinet 32"w x 36"h x 18"d w/1 drawer, 2 doors and 1 shelf. (**Right hand**)



6204—Wall cabinet 25 1/2"w x 30"h x 12"d w/2 doors 6304—Base cabinet 44"w x 36"h x 18"d w/1 drawer & 3 doors 6093—Option: Stainless steel sink w/goose neck faucet

Exam Room Cabinets





6305—Base cabinet 50"w x 36"h x 18"d w/2 drawers, 3 doors and 2 shelves

6096—Option: Work station slide out



6206—Wall cabinet 49 1/2"w x 30"h x 12"d w/3 doors 6306—Base cabinet 50"w x 36"h x 18"d w/2 drawers, 3 doors and 2 shelves.



6096—Work surface pull-out



6207—Wall cabinet 49 1/2" w x 30"h x 12"d w/3 doors

6307—Base cabinet 68"w x 36"h x 18"d w/2 drawers & 4 doors.

6093—Option: Stainless steel sink w/goose neck faucet.

6096—Option: Work station slide out.

6208—Wall cabinet 73 1/2"w x 30"h x 12"d w/5 doors 6308—Base cabinet 74"w x 36"h x 12"d w/3 drawers, 5 doors and 3 shelves

AVAILABLE IN


4 SOLID COLORS &

4 WOODGRAIN COLORS



Walnut











Exam Room Cabinets are very durable and easy to clean. The components are similar to those in our Modular Medical Cabinet line. Modular cabinetry can easily be re-arranged to meet your desired look and efficiency. Important information concerning Options:

6093—Stainless steel sink w/goose neck faucet & strainer basket is available in 24" and 30" base cabinets only. We will make the sink cutout in the counter top when the cabinet is manufactured. To ensure there is no mistake we ask that you submit with your order a photo copy of the group you are ordering and indicate the location of the sink by placing a (S) on the base cabinet pictured. There is no drawer at sink location.

6096—Work surface slide outs are available on 18" & 24" cabinets only. Please place a (W) where desired.

6099—Locks, you may choose to lock 1 or more doors and/or drawers. (note: 2 door cabinets require only 1 lock) Please place an (X) on each door and drawer on your photo copy and submit this information with your order along with the Group number (s) color and please note right or left handed cabinet (s).





Model 6110 - Cabinet w/hutch < (left) **Shown** in Maple w/aluminum pulls This cabinet is built to last and has plenty of room for storage. Makes a great instant office featuring plenty of counter space, storage drawers and cabinets. White counter top only Dimensions: 73"H x 48" W x 19"D

Model 6112 - Storage closet (Right) > Shown in Cherry w/wire pull
This well-built storage closet is perfect for keeping all of your medical materials.
Six shelves and hanging rod included. Two divided vertical spaces make this unit just right to accommodate individual needs.
Dimensions: 71"H x 36"W x 20"D

Locks are available for all cabinets. Order a gang lock for cabinet drawers or single lock which will lock single or double doors.







LOCKERS: Lockable space for cloths and valuables. Each unit includes 2 lockers 36"h x 12"w x 17"d. Orders as many as your space requires then add end panels where needed. You can add additional units anytime.

7600: Locker section 12"w x 17"d x 76"h 7610: End panel 3/4" x 17"d x 74"h

Shown: Maple



7400: Deluxe Exam Table w/adjustable back & footrest 48" to 68"L x 24"W x 32"H Shown: Turquoise / Walnut

7402: Deluxe Exam Table w/adjustable footrest 48" to 68"L x 24"W x 32"H

7404: Deluxe Exam Table w/NO back or footrest 48"L x 24"W x 32"H

















Storage, Lockers, Exam Tables & Carts



6114: Mobile Cabinet w/Hutch Features base cabinet w/casters, 1 door & 2 drawers in Oak finish Hutch w/2 doors and 2 shelves 25"W x 17 1/4"D x 60"H



6116: Treatment Cabinet w/laminated work surface. 2 drawers and 2 door and fixed shelf. Cherry 36"W x 18"D x 36"H



6118: Mobile Utility Cart with 1 drawer, 1 door

Shown: Walnut / wire pull 20"W x 18"D x 30"H



6120: Mobile Utility Cart w/work surface 1 door & 2 drawers. Shown: Oak / wire pulls 25"Wx 17 1/2"D x 32"H



6122: Utility Table or Night Stand w/1 drawer

18"W x 16"D x 24"H



6124: Equipment Cart with 2 drawers.

22"W x 16"D x 36"H (Equipment not included)



6126: Adjustable Cart with 2 adj. shelves.

22"W x 16"D x 36"H



6128: Mobile Utility Cart with 3 drawers

18"W x 19 1/2"D x 30"H



6130: Mobile Cart with 1 door.

Shown: Cherry / alum pull 20"W x 20"D x 31"H



6231: Mobile Cart with one drawer & shelf.

Shown: White / oak pull 20"W x 20"D x 31"H



6232: Mobile Cart with 2 drawers

Shown: Folkstone / oak pulls 20"W x 20"D x 31"H



Modular Cabinetry

PRODUCT DESCRIPTION & TECHNICAL DATA

Modular Medical Cabinetry

Counter top: Made from high-pressure laminates in one-piece construction with the exception to the mitered corner cabinet. Designed with a 3 1/2" back splash and rounded anti-drip front. CARB Phase 2 complaint.

Wood components: Nurse Cabinetry made to compliment Couch products (All wood components) made of 3/4" Melamine with the exception of the drawer sides and cabinet backs which are 1/2". *Funder America Melamine Panels meet or exceed American National Standard A208.1 for 1-M-2 particle board, ALA - NEMA LD3-1995 and meets all requirements of applicable standard CS-236-66 for 1B2 mat-formed particle board. MDF substrates are also used w/melamine panel products. MDF-(medium density fiberboard) and Melamine thermo-fused laminated products are TSCA Title VI complaint.

The above applies to our standard line made from Melamine panels. Should you choose, a full spectrum of Formica or Wilson Art solid colors and wood grains in standard grade laminates are available.

Wood components: Sports Cabinetry made to compliment Table products. Cabinet interiors are made from the same materials as the Nurse Cabinetry, however the exterior components such as the dress panels, doors, and drawer fronts are made from wood veneers, most commonly 7 ply White Maple to insure continuity in grain and color. TSCA Title VI and CARB Phase 2 complaint.

Finish: Available in walnut, cherry, oak or natural finish with a catalyzed water borne top coat for an extremely hard and durable finish. *A gualente* stains and top coats are GREENGUARD certified products from M.L. Campbell wood finishing systems which produce superior resistance to marring, scratching, water, alcohol (200 proof) and household chemicals.

Pulls: Drawer and door pulls are available in oak, maple, aluminum and black or white hooded pulls. (Sports Cabinetry commonly use satin chrome wire pulls which are also available)

Hardware: Hinges are concealed Euro-style Blum series snap-on with 110 degree opening and independent side, depth and height adjustments. Drawer slides are premium quality bottom mount 100 lb. capacity with delron rollers tested to 50,000 cycles. All units have an independent leveling system.

Glass: Glass doors are available on wall cabinets only and are cut from 1/8" acrylic.

Sink: 15" x 15" bar sink with gooseneck faucet, wrist handles and strainer basket. (Not included—supply lines, trap or tail piece.)

Shipping: Depending on size most configurations are shipped on a single sturdy pallet, under one counter top.











Laptop Workstation / Whirlpool Table




Model 9240: without lock >>







Model 9241: with lock >>>> (Laptop computer not included)



This sleek small wall cabinet is typically mounted in the hallways of hospitals, clinics and small rooms where space is limited. Used to securely store and/or charge a Laptop via USB charger (included) and to create a spacious 12" x 18" work surface when opened. Choose our model 9240 without lock or model 9241 which features a Dual Access lock that opens with a key or combination. Use the key to override and/or change the combination which may be shared to limit the need for several keys.

This unit is extremely durable, made from laminated Melamine and features rounded corners and edges.

Available in our standard cabinet colors shown below to compliment our full line of cabinet products.

24"W x 5 3/4"D x 24"H





















7500: Whirlpool Table w/wood frame and laminated top & benches. 48" W x 38" front to back. Shown: Black top w/Oak stand

Available in table finishes only: Cherry / Walnut / Oak or Maple

When ordering please supply the width of your whirlpool tub and floor to rim height.











45" W x 30"D x 30" H
9123: Optional hutch (center)
45" W x 14 1/4" D x 36"H
9124: Optional keyboard pull out
24"W x 16"D

Natural Oak w/aluminum pulls

9122: Desk with drawers on left & right (shown left) 60" W x 30" D x 30" H 9121: Desk with drawers on left (not shown) 45" W x 30" D x 30" H

9120: Desk with drawers on right (not shown) 45" W x 30" D x 30" H

The hutch above the computer desk may also be placed on the 9120 & 9121 desk. (Computer, monitor, etc. not included)



Typical pencil and file drawer used in desk & lateral files.

Pencil drawer - 3/4 extension File drawer - full extension w/folder rail.



9111: Lateral File Cabinet with two drawers 36"W x 18 1/2"D x 31"H

Natural Oak w/hooded pull



9110: Lateral File Cabinet with three drawers 36"W x 18 1/2"D x 43"H

Walnut w/black continuous pull



Bookcase

9100: Book case w/2 shelves 36"W x 13 1/4"D x 30"H Shown: Oak

9101: Book case w/3 shelves 36"W x 13 1/4"D x 43"H

9102: Book case w/4 shelves 36"W x 13 1/4"D x 56"H

Desks are available in Walnut and Oak only.

Bookcases & file cabinets are available in the melamine colors below.

Select Pull style from our sample page 28.



WALNUT.

WALNUT



Office Furinture

PRODUCT DESCRIPTION & TECHNICAL DATA

Office furnishings

Desk top and end panels: Made of 1 1/8" Melamine thermo-fused laminated panels in Oak or Walnut only.

Wood components: All wood components made of 3/4" Melamine with the exception of the drawer sides and cabinet backs, which are 1/2". * *Funder* melamine panels meet or exceed American National Standard A208.1 for 1-M-2 particle board, ALA - NEMA LD3-1995 and meets all requirements of applicable standard CS-236-66 for 1B2 mat-formed particle board. MDF substrates are also used in Melamine panel products. MDF-medium density fiberboard and Melamine thermo-fused laminated product are TSCA Title VI compliant.

Stock colors: white, almond, folkstone, black and fog - Wood grains: oak, maple, walnut and cherry Desks are available in Walnut or Oak finishes only...

Wood components: Sports Cabinetry made to compliment Table products. Cabinet interiors are made from the same materials as the Nurse Cabinetry, however the exterior components such as the dress panels, doors, and drawer fronts are made from wood veneers, most commonly 7 ply White Maple to insure continuity in grain and color.

Finishes: Available in walnut, cherry oak or natural finish with a catalyzed water borne top coat for an extremely hard and durable finish. *A gualente** stains and top coats are GREENGUARD certified products from M.L. Campbell wood finishing systems which produce superior resistance to marring, scratching, water, alcohol (200 proof) and household chemicals. (Pertains to models shown (7300 / 7310 & 7500)

Pulls: Drawer and door pulls are available oak, maple, aluminum and black or white hooded pulls. (Sports Cabinetry Commonly use satin chrome wire pulls which are also available.

Hardware: Hinges are concealed Euro-style Blum series snap-on with 110 degree opening and independent side, depth and height adjustments. Drawer slides are a premium quality bottom mount 100 lb. capacity with delron rollers tested to 50,000 cycles. All units have an independent leveling system.

Shipping: Depending on size most configurations are shipped on a single sturdy pallet.

Customized orders:

Custom quotes are available upon request. Simply submit your millwork plans and/or blue prints.

Thank you for considering WMC for your cabinetry needs.













Space Advantage Table

New Product Introducing our folding Space Advantage Treatment Table

Our *Space A dvantage* wall mount design requires less than 11 cubic feet of space when in a stored position. This product can be used as a Work Station, or with Upholstered Option converts to a Treatment Table. Simple to install and set up, this unit provides a spacious 24 x 72" laminated surface or upholstered top. The table locks securely in the stored position allowing authorized use only. Theses Tables are made of easy to clean durable laminates.

Upholstery option includes: (1) 2" x 24" x 72" snap on padded top, (1) roll of paper w/dispenser & cutter band. This table is also available in a 1" upholstered top version that includes (1) roll of paper w/dispenser & cutter band. The table weights 148 lbs. and will support an approximate weight of 300 lbs.



Work Station Model 7330



only



Work Station w/ upholstery option Model 7340





Upholstered top version



Model 7350

1"

Please specify colors when ordering

Cabinet colors: Cherry, Walnut, Oak, Maple, Almond, White, Folkstone, or Fog. Select upholstery color from samples on page 31.

Vew Produc



9210—Vision & Hearing Screening Cabinet w/acoustic foam and drape. Side table folds down and front table slides into cabinet for storage. Available in White only w/black foam & drape 29"W x 27 5/8"D x 61"H Equipment not included



9360: RX Stool with adjustment ring under seat to adjust height from 19" to 25"

6320



6320— Exam Room Desk Ideal for small spaces Feature storage, drawer and work surface. Shown: Fog w/wire pulls Available in all cabinet colors 21"W x 19"D x 31 1/2"H



4010— Changing Wall Organizer with acrylic dividers 39 1/2" W x 26" H x 13"D Cubical— 12" X 12"

4011 - Changing Cabinet 2 door cabinet with lots of lockable storage Includes pad, restraint and lock 39 1/2" W x 36" H x 22" D Available in all cabinet colors >>>>>





4040– Changing Table

This is an extremely durable cabinet made from select white maple squares and wood veneers. Cabinet features locking storage w/shelf for storing diapers and other baby products. Includes pad, restraint, paper / paper dispenser and lock. For safety, all edges have a rounded profile.

22"W x 52"L x 38"H Available in Table colors: Cherry / Walnut / Oak or Maple only

27



4020– Changing Station with 3—12" x 12" organizers and lockable storage cabinet with shelf.

Includes pad, restraint & lock. 39 1/2"W x 72"H x 22"D Available in all cabinet colors below.









Maple



Almond









WHEN PLACING ORDER PLEASE SELECT MODEL NUMBER - VINYL COLOR CABINET AND/OR TABLE FINISH AND PULL STYLE FROM CHART BELOW Please note: Table finishes are Maple—Oak—Walnut or Cherry only.

Melamine cabinet finishes

| OA | K | CHERRY | MAPLE | WALNUT | ALMOND | FOG | WHITE | BLACK | FOLKSTONE |
|---------------|---|--------|-------|--------|--------|-----|-------|-------|-----------|
| THE THE STATE | | | | | | | | | |
| | | | | | | | | | |
| | | | | | | | | | |
| | | | | | | | | | |

TABLE FINISHES

OAK WOOD VENEER W/SATIN CHROME WIRE PULL

MAPLE WOOD VENEER W/WHITE HOODED PULL

OAK WOOD VENEER W/OAK CONTINUOUS PULL

MAPLE WOOD VENEER W/ALUM. CONTINUOUS PULL

WALNUT WOOD VENEER W/WALNUT CONTINUOUS PULL

OAK WOOD VENEER W/MAPLE CONTINUOUS PULL

CHERRY WOOD VENEER W/CHERRY CONTINUOUS PULL

MAPLE WOOD VENEER W/MAPLE CONTINUOUS PULL

WALNUT WOOD VENEER W/SATIN CHROME WIRE PULL



SAMPLES AVAILABLE BY REQUEST CALL 1-800-757-8055




WOOD FINISH SAMPLES

Walker-McGuire & Companies, Inc. 588 Country Club Rd Pocahontas AR 72455

800 757 8055 email: wmcinc@suddenlinkmail.com

website: wmc-products.com





AGUA SERIES Water Borne Wood Finishes applied to TREATMENT TABLES - TAPING TABLES STATIONS / WHIRLPOOL TABLES & SPORTS CABINETRY.

SEE BACK PANEL FOR PROPER CARE INSTRUCTIONS
A GREENGUARD CERTIFIED PRODUCT FROM:
M. L. CAMPBELL® WOOD FINISHING SYSTEMS



MEDICAL AND SPORTS WOOD FURNISHING ARE FINISHED WITH AGUA SERIES Water Borne Wood Finishes

A Product of M. L. Campbell® Wood Finishing Systems

Superior • Marring • Alcohol (100 proof)

Resistance to: • Scratching • Household chemicals & water

Proper Care Instructions

Remove Dust - Dust buildup from small airborne particles may scratch or dull the surface if not removed correctly. Simply wipe the surface with a cloth dampened with a cleaning polish or mild detergent.

Clean - Fingerprints, cooking fumes, smoking residue, etc., accumulate on any finished surface. None of these contaminants will harm the finish, but they should be periodically removed to restore the finish to its original luster. Just wipe the surface with a cloth dampened with a non-wax containing polish or mild detergent solution. Avoid the use of ammonia based products or silicone oil, as they may cause damage if used over a long period of time.

The finish is durable and resistant to most household spills. However, spills should be wiped up promptly to avoid potential problems.

Excessive prolonged exposure to direct sunlight, high temperatures and high humidity can cause damage to both the finish and the wood itself. These should be avoided.

Following these simple steps will keep your finished piece looking like new for many years.

GREENGUARD Indoor Air Quality Certified®

GREENGUARD Environmental Institute (GEI), in conjunction with the US Environmental Protection Agency (USEPA), has established the most stringent performance-based standards to define products with low chemical and particle emissions for use indoors.

GREENGUARD Certification Program assures the public that VOCs emitted from products like paints and finishes are controlled at the source.

M. L. Campbell finishing products are the first of their kind to achieve GREENGUARD Certification and assure our customers that our products improve indoor air quality and are environmentally conscious.

For GEENGUARD Certified projects, insist on M. L. Campbell's AGUA SERIESTM water borne finishes. Each of the GREENGUARD Certified products is ideal for kitchen, bath or general interior use and delivers the high performance finish you expect from M. L. Campbell.



VINYL COLOR CHART

Walker-McGuire & Companies, Inc. 588 Country Club Rd Pocahontas AR 72455

800 757 8055 email: wmcinc@suddenlinkmail.com website: wmc-products.com









"Naugahyde" Vinyl Cleaning Guide

Cleaning Naugahyde

Advanced BeautyGardTM Protective Finish



Cleaning is a breeze with *Naugahyde* brand fabric and its advanced BeautyGard protective finish. While staining and soiling exposures are common to upholstery fabrics, most stain and soiling may be removed by using the cleaning methods that follow:

For Light soiling:

- 1. A solution of 10% household liquid dish soap with warm water applied with a soft damp cloth will remove most soiling.
- 2. If necessary, a solution of liquid cleanser and water applied with a soft bristle brush. Wipe away the residue with a water dampened cloth.

For Heavier soiling not solved by above method:

1. Dampen a soft white cloth with lighter fluid (naphtha) and rub gently. Rinse with water dampened cloth.

NOTE: Use extreme caution with this method. Complete only in a well-ventilated area away from any open flame.

NOTE: You should try this method on an inconspicuous spot before using it on the original stain/soiling.

For the most difficult stains not removed by the above methods:

- 1. Dampen a soft white cloth with a solution of household bleach (sodium hypochlorite), 10% bleach, 90% water.
- 2. Rub gently.
- 3. Rinse with a water-dampened cloth to remove bleach concentration.
- 4. If necessary, allow a 1:10 diluted bleach solution to puddle on the affected area or apply with a soaked cloth for approximately 30 minutes. Rinse with water-dampened cloth to remove any remaining bleach.

NOTE: You should try this method on an inconspicuous spot before using it on the original stain/soiling.

NOTE: To restore luster, a light coat of spray furniture wax can be used. Apply for 30 seconds and follow with a light buffing using a clean white cloth.

Antimicrobial

Naugahyde with Advanced BeautyGard provides outstanding protection in difficult medical and health care environments. Contains agents effective against bacterial and fungal microorganism.

Please review carefully

CAUTION: Should you desire to use other cleaning methods, carefully try them in an inconspicuous area to determine potential damage to the material. Never use harsh solvents or cleaners which are intended for industrial applications. To clean stained or soiled areas, a soft white cloth is recommended. Avoid use of paper towels.

CAUTION: Cleaning products may be harmful/irritating to your skin, eyes, etc. Use protective gloves and eye protection. Do not inhale or swallow any cleaning product. Protect surrounding area/clothing from exposure. Use in a well ventilated area. Follow all product manufacturer's warnings. Naugahyde cannot be held responsible for damage or injuries resulting from the use or misuse of cleaning products.

ALL WMC UPHOLSTERED PRODUCTS ARE COVERED IN "NAUGAHYDE" A PRODUCT OF UNIFOYAL ENGINEERED PRODUCTS. LLC.





Modular Taping Center design guide

| Project: | | | Order date: | | | | | | | | |
|--------------------------------|---------------------------|---|------------------|--------------------------|-----------|------------|-----------|---------------|--------------|--|--|
| Contact: | | | | Requested delivery date: | | | | | | | |
| Contact phone () | | | Purchase order # | | | | | | | | |
| Uphols nabar Sand / | t ery colors : (ci | ors: (circle one rcle one) Rose / Burgundiew our color op | Orange / Yello | Forest / Dutch | Blue / Re | eg. Blue / | Turquoise | / Jade / Teal | / Chestnut | / Blue Ridge | |
| Budget | | | | | | | | | Сс | st \$ | |
| Product # | | | | | | | | | Use for I | only product numbers lutches & Backs. | |
| CT < Table Extension Product # | | | | | | | | | | Cabinet top Table Extension Product # | |
| Product # | | | | | | | | | | only product numbers Cabinets & Seating | |
| Width | | | | | | | | | То | tal width | |
| Budget | | | | | | | | | Со | st \$ | |

Use the blocks above to draw your design to see how it will look before ordering. **Please note:** Cabinets may join one another however Seating must be placed between two cabinets!

Total project cost \$

Thank you for choosing W M C for your furniture needs.

Please submit this form with your order to insure your design is received and produced as you have illustrated. We will review your order and submit it to production only after we are completely confident that we are all on the same page. If we have questions concerning your order, we will call your contact person at the phone number you have provided to discuss and correct any errors found on this design guide.

W M C



Walker-McGuire and Companies, Inc. 588 Country Club Road Pocahontas, AR 72455